

Interactive Health HT-140 demo video

Written by: Robert Brilliant RB Productions

Client: ICV Digital Media for Interactive Health, Inc.

VIDEO	AUDIO
FADE IN:	(1) FADE UP:
1. Crane shot of model enjoying	Background music theme, in and under
HT 270. Establish and pan into	(1) NARRATOR (VO): Imagine the perfect
set limbo background and	escape from stress for total relaxation in the
dissolve to:	privacy of your own home.
	It's not just a fantasy anymore
2. Crane beauty shot of empty	(2) Meet your personal masseuse, the HT-
HT-140 in classic limbo set.	140 robotic massage chair, from Interactive
Animate title over shot: <i>HT-140</i>	Health, the world leader in massage technology.
from Interactive Health.	
3. MCU back of chair. Key	(3) The HT-140 is built upon patented
animated HTT logo. Key insert	Human Touch Technology that closely replicates
roller assembly in motion.	therapeutic techniques used by chiropractors and
	massage professionals.
4. Montage: Top view looking	(4) This advanced technology employs a

down on massage motor	three-motor mechanism that follows the spine's
mechanism in motion.	natural curve.
	To this soothing motion we've added a
Montage showing soothing	unique quad roller assembly that performs like
motion. Key: Performs like an	an expert masseuse.
expert masseuse.	
5. Montage: MS and MCU of	(5) The result is a multi-dimensional robotic
model enjoying massage in HT-	massage that feels remarkably effective, relaxing,
140. Motion shot.	and human!
6. MS: another model in chair.	(6) Only chairs with Human Touch
Key: All Major Massage Modes	Technology, like the HT-140, can replicate all
	major massage modes. These include:
CU rolling motion, Key: •	-Rolling, which gently relieves tension and
Rolling	loosens muscles along the spine,
Split screen showing hands-on	-Kneading, a circular, pressing motion that
rolling massage. Cut to:	goes deeper to lift and stretch muscles for
Split screen: hands-on kneading	improved circulation,
and roller kneading with title: •	-Percussion, a rapidly tapping motion that
kneading, cut to: front angle CU	flexes spinal joints, relieves pressure, and
of kneading rollers. Cut to:	invigorates the entire back, and
Split screen: hands-on percussion	-Compression, a mode exclusive to Human
and roller percussion with title: •	Touch Technology, that firmly pushes muscles on
percussion, cut to: front angle	each side of the spine, for improved mobility and
	posture.

	1
CONT. 6 CU of percussion rollers. Cut to: Split screen: hands-on compression and roller compression with title: • compression, cut to: front angle CU of compression rollers.	
7. MS model in chair. Move to	(7) The HT-140 can even combine kneading
MCU. Key: Kneading and	and percussion for the most intense, deep-tissue
Percussion	massage you've ever experienced in a robotic
	chair.
8. MCU: foot and calf massager	(8) The HT-140 also comes with this multi-
in use. Key: Foot & Calf	speed foot and calf massager that delivers
Massager	reflexology-style benefits to the entire body.
9. CU: neck massager in action,	(9) A penetrating neck massager completes
model lays back into chair, fast	the HT-140's head-to-toe massage system.
zoom out to MWS.	
10. CU, controller in person's	(10) This hand-held controller is your HT-
hand.	140 command center.

	·
11. ECU side buttons. Highlight as mentioned.	(11) With these buttons on the controller's side, you can adjust shoulder height, and select one of three fifteen-minute auto massages for full, upper, or lower back.
12. Composited shot: MS model in chair, ECU width adjust button, with title: Customizable Width	(12) The width button adjusts the massage roller width for customizable comfort.
13. Composited shot: MWS person receiving stretch massage. ECU stretch button, with title: Stretch Mode	(13) The stretch button activates our unique stretch mode, which emulates traction by gently stretching and uncompressing the spine.
14. Composited shot: MWS person in chair. ECU restore button, with title: <i>Restore</i> . Lose title and show restore mode in action.	(14) When you're through massaging, the restore button brings the chair to a full upright position and stows the footrest.

15. ECU stop button. Key: Stop	(15) To disengage any chair function at any time, simply press the stop button.
16. MS Chair on set, with model reading a book or listening to iPod.	(16) Even when not massaging, the HT-140 provides remarkable comfort and relaxation.
17. WS Chair reclining to 170 degrees, insert box showing remote controller. Key text: 170 Degree Recline	(17) Through the controller, the HT-140 can be reclined 170 degrees.
18. Split screen: CU removable softening pad, head/neck pillow	(18) For added comfort, there's a removable softening pad for the back and an extra thick head and neck rest pillow.
19. Montage: Build of HT-140 stock photos in home environment.	(19) More than just a great massage chair, the HT-140 is beautiful furniture you'll be proud to own. Its exquisite design is quality crafted with sumptuous leather and faux leather in popular black and cashew colors.
20. Montage: CUs oak armrests and base.	(20) Rich oak wood base and armrests complete this elegant fit and finish.

21. Model sitting in chair. Crane shot. Key: Sitting is Believing	(21) At Interactive Health, we firmly believe that "sitting is believing." Spend just a few minutes in the HT-140 and you'll be convinced: robotic massaging at long last has "the human touch."
21. Beauty shot of chair with model, go to silhouette lighting. FADE TO BLACK	(21) The HT-140 from Interactive Health, another quality product from the world leader in massage technology. Music up and Fade Out
END	